



INTERVENTIONAL  
VASCULAR  
DIAGNOSTICS  
AND THERAPY

# SeQuent<sup>®</sup> Please NEO

CLINICALLY PROVEN DRUG COATED BALLOON CATHETER

# SeQuent<sup>®</sup> Please NEO

CUTTING-EDGE DRUG COATED BALLOON CATHETER


## THE SECOND GENERATION DCB

### Outstanding performance:

- Advanced crossing properties
- Improved pushability
- Hydrophilic shaft coating
- Reduced balloon wall thickness

### Clinically proven indications:

- In-stent restenosis
- De novo
- Small vessel disease
- Bifurcations




**OVER 28 CLINICAL TRIALS WITH OVER 5,900 ENROLLED PATIENTS**

## IMPLANT-FREE WITH SeQuent® Please NEO

No stent-related complications and only **1-month DAPT** for the treatment with DCB-only


### Clinically proven Paclitaxel and Iopromide coating

The matrix coating of Paclitaxel and Iopromide ensures the effective drug release into the vessel wall.



### Homogenous drug delivery<sup>1-4</sup>

Only a "single shot" drug delivery with SeQuent® Please NEO is needed to ensure a sustained anti-proliferative effect. A short inflation time of only 30 seconds proved to be sufficient to inhibit cell proliferation.<sup>2</sup>



Stent struts of a DES lead to an inhomogeneous drug distribution pattern. About 85% of the vascular wall is not covered by the struts resulting in low drug tissue level.<sup>1</sup>




Homogenous drug distribution with SeQuent® Please NEO.<sup>3</sup>

## PROVEN LATE LUMEN ENLARGEMENT

SeQuent® Please NEO supports the inherent mechanism of natural vessel restoration and leads to late lumen enlargement

### Clinical trial to study late lumen enlargement of de novo lesions after DCB-only<sup>5</sup>



Angiographic Measure	Minimal Lumen Diameter in mm
Pre-treatment	0.81 ± 0.47
Post-treatment	1.75 ± 0.58
4-month follow-up	1.91 ± 0.55
p-value pre vs. post	< 0.001
p-value post vs. 4-month follow-up	< 0.001

Late lumen enlargement after 4 months

**+ 0.16 mm**

<sup>1-4</sup> Axel DJ et al. Circulation 1997; 96: 636-45. | Hwang CW et al. Circulation 2004; 104: 600-5. | Scheller B et al. Circulation 2004; 110: 810-4. | Scheller B et al. Heart 2007; 93: 539-41.

<sup>5</sup> Kleber F et al. Clin Res Cardiol 2015; 104: 217-25.

# SeQuent<sup>®</sup> Please NEO

DCB-ONLY TREATMENT

## ADVANTAGES OF DCB-ONLY

### No unnecessary stent implantation


- No inflammation due to a foreign body implant
- No risk of stent thrombosis
- No stent-related limitations for further treatment
- No stent edge effect

### Efficacy of DCB

- Enable positive remodeling
- Keep natural vessel vasomotion
- Only 1-month DAPT: Cost efficacy studies ongoing

DCB-only provides the standard of care for all patients with high bleeding risks and atrial fibrillation<sup>⑥</sup>

## METHODOLOGY<sup>⑦</sup>




⑥ Valgimigli M et al. European Heart Journal 2018; 39(3): 213-60.

⑦ Kleber F et al. Clin Res Cardiol 2013; 102: 785-97.

# GO IMPLANT-FREE

**BASKET-SMALL 2: Randomized clinical trial for DCB-only vs. DES in de novo lesions (small vessel disease)<sup>8</sup>**




**Primary endpoint: MACE at 12-month follow-up in %**

DES (Xience/ Taxus® Element™)	7.54
DCB-only (SeQuent® Please NEO)	7.57
p-value	0.92

**DCB-only is non-inferior to DES in de novo lesions up to 3 mm**

**OCTOPUS II: Clinical trial using OCT to evaluate the use of DCB without stenting in de novo lesions<sup>9</sup>**



**Primary endpoint: Late Lumen Loss at 6-month follow-up in mm**

DES (Xience) <sup>10</sup>	0.16 ± 0.15
DCB-only (SeQuent® Please)	-0.13 ± 0.44
p-value	< 0.05

**DCB-only achieves long-term late lumen gain contrary to DES**

<sup>8</sup> Jeger R et al. The Lancet 2018; 392(10150): 849-56.

<sup>9</sup> Poerner T et al. Clin Res Cardiol 2017; 106: 18-27.

<sup>10</sup> Poerner T et al. CCI 2014; 7(6): 760-7.

# SeQuent® Please NEO


CLINICALLY PROVEN INDICATIONS

## IN-STENT RESTENOSIS


**Patient:** Male, 55 years

**Indication:** ISR of BMS (3.5 x 15 mm) implanted 2 years ago


**Procedure:** Pre-dilation 3.5 x 15 mm PTCA balloon  
DCB-only SeQuent® Please (3.5 x 20 mm) proximal lesion  
DCB-only SeQuent® Please (3.5 x 15 mm) distal lesion



Pre-treatment



Post-treatment



4-month follow-up

Drug coated balloons are recommended for the treatment of in-stent restenosis (BMS or DES) by the ESC Guidelines<sup>①</sup>

I

A

507-511,524

## DE NOVO LESION

**Patient:** Female, 67 years

**Indication:** De novo stenosis of obtuse marginal branch


**Procedure:** Pre-dilation 2.5 x 15 mm PTCA balloon  
DCB-only SeQuent® Please (2.5 x 20 mm)



Pre-treatment



Post-treatment




4-month follow-up

## BIFURCATION


**Patient:** Male, 54 years

**Indication:** Stenoses of mid circumflex artery (CX) and its posterolateral branch (PL-CX)


**Procedure:** Pre-dilation 2.5 x 20 mm PTCA balloon of CX  
DCB-only SeQuent® Please (3.0 x 15 mm) of PL-CX  
DCB-only SeQuent® Please (3.0 x 20 mm) of CX



Pre-treatment



Post-treatment




4-month follow-up

<sup>①</sup> Windecker S et al. European Heart Journal 2014; 35: 2541-619.

Balloon Diameter	Balloon Length							Nominal Pressure	Rated Burst Pressure
	10 mm	15 mm	20 mm	25 mm	30 mm	35 mm	40 mm		
2.0 mm	5023200	5023210	5023220	5023230	5023240	5023250	5023260	6 atm	14 atm
2.25 mm	5023201	5023211	5023221	5023231	5023241	5023251	5023261	6 atm	14 atm
2.5 mm	5023202	5023212	5023222	5023232	5023242	5023252	5023262	6 atm	14 atm
2.75 mm	5023203	5023213	5023223	5023233	5023243	5023253	5023263	6 atm	14 atm
3.0 mm	5023204	5023214	5023224	5023234	5023244	5023254	5023264	6 atm	14 atm
3.5 mm	5023206	5023216	5023226	5023236	5023246	5023256	5023266	6 atm	14 atm
4.0 mm	5023207	5023217	5023227	5023237	5023247	5023257	5023267	6 atm	14 atm

Technical Data	
Proximal shaft	1.9 F
Distal shaft	2.5 F
Usable length	145 cm
Balloon crossing profile	0.033" - 0.037"
Lesion entry profile	0.016"
Guiding catheter compatibility	5 F standard guiding catheter
Guidewire compatibility	0.014"
Rated burst pressure [RBP]	14 atm
Nominal pressure [NP]	6 atm





**Distributor**

B. Braun Melsungen AG | Vascular Systems | Sieversufer 8 | 12359 Berlin | Germany  
 Phone +49 30 568207-300 | Fax +49 30 568207-210 | [www.bbraun.com](http://www.bbraun.com)

Manufacturer acc. to MDD 93/42/EEC of SeQuent® Please NEO is the B. Braun Melsungen AG, Carl-Braun-Str. 1, 34212 Melsungen, Germany.

The product trademark 'SeQuent' and 'Coroflex' are registered trademarks of B. Braun Melsungen AG. 'NSE Alpha' is a trademark of Goodman Co. Ltd., Japan. 'Taxus' is a registered trademark of Boston Scientific Corporation.

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.